



Investigation of Discrepant ABO Blood Grouping Results from an Autoanalyzer

Young Ae Lim , M.D., Ph.D.^{1,2}, Seo-Jin Park , M.D., Ph.D.^{1,2}, and Hyun Soo Cho , M.S.²

¹Department of Laboratory Medicine, Ajou University School of Medicine, Suwon, Korea; ²Department of Laboratory Medicine, Ajou University Hospital, Suwon, Korea

Background: A paucity of studies evaluating failed cases of ABO grouping using autoanalyzers exists. We investigated autoanalyzer rejected cases, including serologically suspicious ABO subgroups and discrepant ABO blood grouping results from Erytra Eflexis (Grifols, Spain), to demonstrate efficient use of autoanalyzers for ABO grouping.

Methods: Samples requested for ABO grouping throughout 2020 were tested using two Eflexis instruments and standard ABO RhD and reverse grouping cards. Neonatal cards were not used. When necessary, a conventional tube technique (TUBE) was used to resolve rejected/discrepant Eflexis ABO grouping results.

Results: The overall sample rejection rate (RR) was 3.2% (628/19,466), 1.3% of which were due to various error flags and 1.9% for discrepant results. Cases from neonates ≤ 1 year old accounted for 35.3% of the rejected cases based on Eflexis results. The ABO groups with the highest and lowest RR (excluding neonates) were A and O, respectively. The 628 samples resulted in 682 rejections, which were frequently associated with reverse grouping, including 28.4% against A₁ and 54.5% for B red cells. Among 14 serologically weakened A and/or B blood groups, six A₂B_w and two AB_w, which had been missed by Eflexis, were detected using TUBE and our follow-up laboratory criteria.

Conclusions: The ABO group and a proportion of neonatal samples influenced the RR due to weak reverse grouping reactivity, especially toward B red cells. Confirmatory ABO grouping by TUBE in a new patient and/or extra rejection criteria for forward grouping are needed to detect cis-AB, which is relatively common in Korea.

Key Words: ABO blood group system, Autoanalyzer, Instrument, Rejection

Received: January 2, 2022

Revision received: March 4, 2022

Accepted: June 15, 2022

Corresponding author:

Young Ae Lim, M.D., Ph.D.
Department of Laboratory Medicine, Ajou University School of Medicine, 206 World cup-ro, Yeongtong-gu, Suwon 16499, Korea

Tel: +82-31-219-5786

Fax: +82-31-219-5778

E-mail: limyoung@ajou.ac.kr



© Korean Society for Laboratory Medicine

This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (<https://creativecommons.org/licenses/by-nc/4.0>) which permits unrestricted non-commercial use, distribution, and reproduction in any medium, provided the original work is properly cited.

INTRODUCTION

Blood bank autoanalyzers can save technician time and prevent various errors associated with manual testing. Since August 2015, the Korea Health Insurance Review & Assessment Service instituted a higher reimbursement fee for automated ABO blood grouping than for manual tests. Since this change, the number of hospitals that routinely use autoanalyzers for ABO blood grouping has significantly increased.

When results from blood bank autoanalyzers are rejected be-

cause of failure to determine ABO blood groups due to discrepant results or serologically suspicious ABO subgroups, the laboratory must have a manual process in place to investigate the problem. Thus, although blood bank autoanalyzers may have several advantages, the constant need for follow-up investigations may increase the laboratory workload and risk of errors. Therefore, it is necessary to identify the characteristics of rejected cases and develop a strategy to minimize their occurrence.

Erytra Eflexis (Eflexis; Grifols, Barcelona, Spain) is an immu-

nohematology autoanalyzer suitable for pre-transfusion routine testing in small-to-medium sized laboratories [1, 2]. Ajou University Hospital, Suwon, Korea, was the first to introduce Eflexis analyzers in Korea in 2019. There are several reports on the performance and reliability of Erytra or Erytra Eflexis in pre-transfusion testing [1-4]. However, there is a paucity of comprehensive studies evaluating failed cases of ABO grouping using Eflexis or autoanalyzers from other companies [5-9].

Cis-AB, a rare ABO subgroup, is caused by a gene mutation resulting in a single glycosyltransferase enzyme with dual A and B glycosyltransferase activities [10]. The typical phenotype of *cis*-AB is A₂B₃. The *cis*-AB phenotype is the most frequent ABO subgroup in Korea (0.0354%) [11]. It also occurs more frequently in Northeast Asia than in other regions [12]. During our evaluation of Eflexis during 2019, we found several cases where the analyzers had misidentified the A₂B_w phenotype (presumed to be *cis*-AB) as normal AB groups. Similar results of ABO mistyping of the *cis*-AB group by other autoanalyzers have been reported [13]. Subsequently in 2019, our laboratory requested that Grifols make modifications to the autoanalyzer, such that it is able to more clearly differentiate 4+ and 3+ hemagglutination reaction scores when using anti-B anti-sera, to allow for improved differentiation of the A₂B_w phenotype from that of the normal AB group. Consequently, the manufacturer modified the Eflexis software, and an upgraded version (v.1.2.2) was released in August 2020.

Since the introduction of the Eflexis analyzers at Ajou University Hospital, we have included extra criteria to prevent automatic interfacing of discrepant ABO grouping results to the Laboratory Information System (LIS) and to carry out further investigations using the conventional tube technique (TUBE). Further review was required to manage the accumulated data from the additional rejection criteria to detect serologically suspicious ABO subgroups, including the *cis*-AB phenotype.

The aim of this study was to inform laboratories on how to use immunohematology autoanalyzers more efficiently for ABO grouping by reviewing cases with a serologically weakened A or B antigen to detect ABO subgroups or analyzing rejected cases of ABO grouping that were not automatically interfaced by the Eflexis analyzer.

MATERIALS AND METHODS

Subjects

This retrospective study was approved by the Institutional Review Board of Ajou University Hospital (AJIRB-BMR-SMP-20-299).

Samples from all patients who requested ABO blood grouping between January 1 and December 31, 2020, were tested using either of the two Eflexis analyzers. EDTA-anticoagulated whole blood samples were used for ABO grouping. ABO blood grouping results, patient's ward, test request time, ABO grouping hemagglutination reaction strengths, or any flag results for forward or reverse ABO grouping were retrospectively extracted from both Eflexis analyzers as an Excel file (Microsoft 2019, Microsoft, Redmond, WA, USA). For the rejected cases, patient age and sex were additionally extracted. Babies ≤12 months of age were classified as neonates.

To evaluate the ability of the Eflexis analyzers to detect serologically suspicious ABO subgroups, samples with weakened A and/or B antigen reactions identified in our laboratory between January 2020 and April 2021 were also included in the study. These samples were confirmed by TUBE, not by molecular testing. Aw was defined as serologically weakened group A other than A₁ and A₂, and Bw was defined as serologically weakened group B.

ABO grouping using the Erytra Eflexis analyzer

ABO grouping was performed using the Erytra Eflexis instrument and DG Gel 8 ABO/Rh cards (Diagnostic Grifols, Barcelona, Spain), according to the manufacturer's recommendations. Neonatal cards were not used. The ABO cards consisted of eight columns, *i.e.*, anti-A, anti-B, anti-AB, anti-D, anti-D^{VI} clone, control, and two empty wells for reverse grouping using 0.8% A₁ and B reagent red blood cells (RBCs) (Serigrup Diana A₁/B, Diagnostic Grifols), respectively. We used 0.8% patient RBCs suspended in low ionic salt solution for ABO forward grouping and RhD determination. Automated hemagglutination reaction strengths were graded as –, 1+, 2+, 3+, or 4+. Alert flags, such as “trace,” were returned by the Eflexis analyzer when the reaction could not be graded or interpreted. For comparison of the reaction strength between the autoanalyzers and TUBE, “–” and “trace” reactions were considered 0 and 0.5, respectively.

Each sample was classified into a specific ABO group according to the results of forward and reverse grouping interpreted by the system's software, and the final ABO grouping result was automatically interfaced to the LIS. When the autoanalyzers could not classify the ABO group, the results were withheld by the analyzer and not automatically interfaced to the LIS. Such results were identified as “Eflexis rejections.”

ABO blood group confirmation for rejected cases from the analyzers or serologically suspicious ABO subgroups was determined in the following four ways. 1) For patients who previously

underwent ABO grouping in our hospital, the historical ABO group was considered the final ABO group, except for patients who had undergone an ABO-mismatched peripheral blood stem cell transplantation (PBSCT). 2) For new patients, forward and reverse ABO grouping using TUBE was performed, except for neonates, for whom only forward grouping was performed. 3) For cases with discrepant results from previous ABO blood grouping results or discrepant results between forward and reverse grouping, or a suspected ABO subgroup, the ABO group was determined by repeating the ABO blood grouping using the autoanalyzers and TUBE, as well as by reviewing the patient's diagnosis. When necessary, additional tests, including antibody screening and identification, were performed using the autoanalyzers. Further tests using manual methods were performed, including the saline replacement method, washing RBCs with warm saline, and testing RBCs with anti-A₁ (Lorne Laboratories Limited, Lower Earley, UK) and/or anti-H lectins (Lorne Laboratories Limited) using TUBE. 4) If it was difficult to determine the ABO group even after TUBE and additional tests, the result was reported as "Unclassified ABO group."

ABO blood grouping using TUBE

ABO blood grouping using TUBE was performed following the laboratory's standard operating procedure and manufacturer's recommendations for all new patients, irrespective of the autoanalyzer results. IgM monoclonal anti-A and anti-B anti-sera (Shinyang Chemical, Seoul, Korea; raw materials from Bioscot, Millipore, Livingston, UK) were used for forward grouping. For reverse grouping, pooled A₁ and pooled B RBCs (7–10 random donors each for A₁ and B RBCs) suspended in 2%–5% saline were used. Incubation at room temperature for 15 minutes was standard for TUBE-based reverse grouping. The hemagglutination reactions for forward and reverse TUBE-based grouping were scored as follows: 0: negative, 0.5+: weak, 1+, 2+, 3+, or 4+ [14].

Study design

Rejected, non-interfaced cases from the autoanalyzers were classified into the following two groups. 1) Eflexis rejection group: failed analyzer ABO blood grouping and failed automatic interfacing because of the presence of any flag, such as "trace," "double population (DP)," or "?" (defined as the Flag group) or a discrepancy between forward and reverse grouping without any flag (defined as the Discrepant result group). 2) Laboratory rejection group: extra criteria of our laboratory, including cases with ABO hemagglutination reaction scores from 1+ to 3+ for

forward grouping or 1+ for reverse grouping, preventing automatic interfacing, although an ABO group was determined by the analyzers (Fig. 1). The extra criteria were used for the detection of A₂B_w with weak reaction of anti-B antibodies, which was typical for the *cis*-AB group.

The rejection rate (RR) for non-interfaced cases was determined as the percentage of the number of rejected cases divided by the total number of automatically interfaced cases plus rejected cases from the autoanalyzers during 2020. Cases with identical results obtained from the same sample after retesting within 10 minutes were excluded from the rejected cases and were defined as "same repeated cases." The frequency of rejected results was defined as the percentage of the number of rejected results divided by the total number of cases with rejected results from the autoanalyzers during 2020.

Statistical analysis

The chi-square test was used to compare sample RRs among ABO blood groups. The paired *t*-test was used to compare the average strength of reverse grouping reactions between TUBE and the autoanalyzers. Statistical analysis was performed using Excel (Microsoft), and $P < 0.05$ was considered statistically significant.




Fig. 1. Numbers of automatically interfaced and non-interfaced cases. The non-interfaced cases were divided into various flags (Flags), discrepant cell and reverse grouping (Discrepant results), or extra arbitrary laboratory criteria (Laboratory rejections) from two Eflexis analyzers. The rejection groups were further divided into ≤ 1 year to describe neonates versus patients older than 1 year old (> 1 year). *Numbers of cases excluding those with 37 discrepant results between forward and reverse ABO grouping due to previous ABO-mismatched peripheral blood stem cell transplantation.

RESULTS

The sample RR was 3.2% (628/19,466) in the Eflexis rejection group and 1.7% (331/19,466) in the Laboratory rejection group (Fig. 1). The Eflexis rejection group comprised 249 Flag group

cases (1.3%) and 379 Discrepant result group cases (1.9%). Thirty-seven cases were excluded from ABO grouping due to previous mismatched PBSCT, causing discrepant results between forward and reverse grouping. The remaining 591 cases in the Eflexis rejection group could be classified into an ABO

Table 1. Interfaced and non-interfaced cases according to the ABO blood group during 2020

Cases	Age	ABO blood group, n (%*)					Total
		A	B	O	AB		
Interfaced	All	6,176	5,076	5,065	2,190		18,507
Eflexis rejections	≤ 1 yr [†]	104 (1.6)	50 (0.9)	62 (1.2)	6 (0.3)		222 (1.1)
	> 1 yr [†]	218 (3.3)	112 (2.1)	18 (0.3)	21 (0.9)		369 (1.9)
	All [†]	322 (4.8)	162 (3.0)	80 (1.5)	27 (1.2)		591 (3.0)
Laboratory rejections [†]	All	185 (2.8)	88 (1.7)	19 (0.4)	39 (1.7)		331 (1.7)
Total rejections [†]	All	509 (7.6)	250 (4.7)	99 (1.9)	66 (2.9)		922 (4.7)

Eflexis rejections include non-interfaced cases with failed ABO blood grouping on the Eflexis analyzers because of the presence of various flags or discrepancies between forward and reverse grouping without any flags. Laboratory rejections include cases acquired from the extra arbitrary criteria applied in our laboratory, including 1+ to 3+ hemagglutination reaction strength for forward ABO grouping or 1+ for reverse ABO grouping. This prevented automatic interfacing of results, although the ABO blood group was determined by the autoanalyzers.

*% indicates rejection rates; [†]P<0.001 using chi-square test.

Table 2. Eflexis rejection results in 628 cases according to the ABO blood group and ABO grouping reagents

ABO Grouping reagent	ABO group	No. of cases	Hemagglutination reaction strength					Flag			Rejected results	
			4+	3+	2+	1+	0	DP	trace	?	N	%
Anti-A	A	322	290	0	0	0	0	32	0	0	32	4.7
	B	162	0	1*	0	0	158	3	0	0	4	0.6
	O	80	0	0	0	0	80	0	0	0	0	0
	AB	27	14	0	0	0	0	13	0	0	13	1.9
	NA	37	0	0	0	0	0	14	0	0	14	2.1
Anti-B	A	322	0	0	0	0	320	0	2*	0	2	0.3
	B	162	131	7*	0	0	0	24	0	0	31	4.5
	O	80	0	0	0	0	80	0	0	0	0	0
	AB	27	21	0	0	0	0	6	0	0	6	0.9
	NA	37	0	0	0	0	0	14	0	0	14	2.1
A _i RBCs	A	322	0	0	2*	2*	312	0	5*	1	10	1.5
	B	162	14	16	5	2	100 [†]	0	22	3	125	18.3
	O	80	7	10	2	11	40 [†]	0	9	1	50	7.3
	AB	27	0	1*	1*	2*	21	0	1*	1	6	0.9
	NA	37	0	0	0	0	0	0	3	0	3	0.4
B RBCs	A	322	15	11	3	7	197 [†]	0	76	13	286	41.9
	B	162	0	3*	1*	3*	149	0	5*	1	12	1.8
	O	80	3	9	2	2	53 [†]	0	10	1	64	9.4
	AB	27	0	4*	0	3*	17	0	2*	1	10	1.5
	NA	37	0	0	0	0	0	0	0	0	0	0

*Unexpected atypical reactions of ABO grouping; [†]Lack of expected reverse reactions of ABO grouping.

Abbreviations: DP, double population; NA, ABO blood grouping result not available because of ABO-mismatched PBSC transplantation.

group. Neonates accounted for 35.4% (222/628) of the Eflexis rejected cases and 270/959 (28.1%) of the total rejected cases.

The sample RR significantly differed according to the ABO group in the Eflexis rejections ($P < 0.001$), laboratory rejections ($P < 0.001$), and total rejections ($P < 0.001$). Blood group A showed the highest RR, and group O the lowest. Group AB had the lowest RR in neonates (0.3%) (Table 1).

The results (N=682) of the Eflexis rejection group from 628 rejected samples are shown in Table 2. According to the ABO blood grouping reaction, 682 rejected results resulted from anti-A (9.2%) and anti-B anti-sera (7.8%), and reactions against A₁ (28.4%) and B RBCs (54.5%). Rejections for ABO blood group-

ing by the autoanalyzers were caused by a lack of the expected reverse reaction or an unexpected reaction associated with cold auto-antibodies, allo-antibodies, or rouleaux formation (Table 2).

There were 263 flags (106 “DP,” 135 “trace,” and 22 “?” flags) among the 249 flagged cases (Fig. 1, Table 2). “DP” flags were reported only in forward grouping, whereas “trace” and “?” flags were also reported in reverse grouping, except for two cases of “trace” in forward grouping. While “DP” occurred in 26.4% (28/106) of results following ABO-mismatched PBSCT, this flag was also observed in cases with serologically weakened A and/or B groups, neonates, and patients who received emergency transfusions of O RBCs. The two cases with a “trace” flag in forward

Table 3. Reverse grouping reaction strengths for three methods in 112 cases of new patients

Comparison parameter	B RBCs (N=81*)			A RBCs (N=35*)		
	Eflexis	IS TUBE	15-minute TUBE	Eflexis	IS TUBE	15-minute TUBE
Mean ± SD	0.64 ± 0.75	0.99 ± 0.50	1.88 ± 0.56	0.86 ± 1.23	1.11 ± 1.04	1.94 ± 0.94
P vs. Eflexis		<0.001	<0.001		0.0049	<0.001
P vs. IS			<0.001			<0.001

*Four blood group O samples are included.

Abbreviations: SD, standard deviation; TUBE, conventional tube technique; IS, immediate spin; 15-minute, 15-minute incubation at room temperature.




Fig. 2. Examples of different hemagglutination reaction strength interpretations by Eflexis using anti-B anti-sera before (A, case 6 in Table 4) and after (B, case 10 in Table 4) the software upgrade to v.1.2.2. Starting from the left side, the first three columns are for forward grouping, the next two columns are for D typing, then a column for forward grouping control, and the last two columns are for reverse grouping. Similar reaction strengths using anti-B anti-sera were observed before and after the software upgrade. However, while the analyzers had previously determined the anti-B reaction strength as 4+ with software v1.2.1 (red-outlined rectangle, A), it was graded as a 3+ reaction following the software upgrade to v1.2.2 (red-outlined rectangle, B). Examples of expected 4+ reaction strengths (A and B) toward anti-A anti-sera are also shown.

Table 4. Hemagglutination reaction strengths in Eflexis and TUBE for ABO grouping in 14 patients with serologically weakened A or B antigens

Case number	Sex	Age	Eflexis			IS TUBE				15-minute TUBE					
			ABO	anti-A	anti-B	A cell	B cell	ABO	anti-A	anti-B	A cell	B cell			
1	F	57	NA	DP	0	0	4	Aw*	1	0	0	4	0	4	
2	M	65	NA	0	DP	4	0	Bw	0	2	4	0	ND	4	0
3	F	66	NA	0	DP	3	0	Bw	0	1	2	0	ND	3	0
4	M	73	NA	0	DP	4	0	Bw	0	3	4	0	ND	4	0
5	F	57	NA	0	DP	4	0	Bw	0	2	4	0	ND	4	ND
6	M	26	AB	4	4	0	0	A ₂ Bw	4	2	0	0	0	4	0
7	M	59	AB	4	4	0	0	A ₂ Bw	4	2	0	0	0	4	1
8	F	67	AB	4	3	0	0	A ₂ Bw	4	2	0	0	0	4	0
9	M	49	AB	4	3	0	0	A ₂ Bw	4	2	0	0	0	4	1
10	F	85	AB	4	3	0	0	A₂Bw	4	2	0	0.5	0	4	1
11	M	56	AB	4	3	0	0	A₂Bw	4	2	0	0	0	3	0.5
12	M	59	AB	4	1 [†]	0	0	ABw	4	Weak	0	0	4	1	0
13	F	43	AB	4	2	0	0	ABw	4	Weak	0	0	4	Weak	0
14	F	37	NA	4	DP	0	0	ABw	4	1	0	0	4	Weak	0

The eight ABO blood groups highlighted in bold were determined following the software upgrade (v1.2.2) of the Eflexis analyzers.

*Case 1 was a patient diagnosed as having AML whose blood group returned to normal blood group A after three months; †A case with initially 1+ strength but showed a negative result in a repeat test using Eflexis.

Abbreviations: TUBE, conventional tube technique; IS, immediate spin; 15-minute, 15-minute incubation at room temperature; ABO, ABO blood group; DP, double population; NA, not available; ND, not determined; Weak, weak reaction strength < 1+.

grouping were group A and were considered nonspecific reactions in the forward group using anti-B anti-sera. The “trace” and “?” flags represented weak reactions in reverse grouping, and none of these were false reactions. Forty-nine percent (129/263) of the flags were observed in group A. The frequencies of flagged cases were 16.7% (37/222) for patients aged ≤ 1 year old and 55.7% (226/406) for those > 1 year old.

In the Laboratory rejection group, there were 336 rejected results from 331 patients. The rejected results comprised 55 results with a 3+ reaction strength, one result with a 1+ reaction strength in forward grouping, and 280 results with a 1+ reaction strength in reverse grouping, including 190 results from reverse grouping using B RBCs. One hundred eighty-six (55.4%) results were from group A, 88 (26.2%) from group B, 39 (11.6%) from group AB, and only 19 (6.8%) from group O. Except for one case, group AB rejections were because of reaction strengths $\leq 3+$ for forward grouping using anti-B anti-sera.

Reverse grouping using TUBE was performed in 112 patients. The average hemagglutination reaction strengths for reverse grouping by TUBE with immediate spin and following a 15-minute incubation at room temperature were significantly higher than those for the autoanalyzers (Table 3).

During the 16-month study period, 14 cases of serologically suspected ABO subgroups were observed: one A_w , four B_w , three AB_w , and six A_2B_w (Table 4). Thirteen cases were new patients who required retesting for ABO blood grouping using TUBE. Case No. 1 was a patient with AML whose ABO grouping returned to normal blood group A after three months. Because of a DP flag in forward grouping, the results were not automatically interfaced for six cases of A_w , B_w , and AB_w , which were easily detected. However, six cases of A_2B_w and two cases of AB_w were determined as group AB and would have been automatically interfaced by the autoanalyzers had our laboratory rejection criteria not been in place. The false but strong reactions in these cases were detected by laboratory staff by performing repeat TUBE not only for new patients but also for samples with reaction strengths of 1+ to 3+. All six A_2B_w cases showed 2+ reaction strengths and mixed-field reactions in TUBE-based forward grouping using anti-B anti-sera, which are typical results for *cis*-AB. Anti-B antibody was detected after a 15-minute incubation at room temperature in four cases of A_2B_w , whereas the autoanalyzers did not detect any anti-B antibody. Using the Eflexis software v1.2.1, two out of four cases (case Nos. 6 and 7, Table 4) showed 4+ hemagglutination in forward grouping using anti-B anti-sera, although the reaction strength was weaker than the typical 4+ strength for a normal AB group. Following the software upgrade

(v1.2.2 or later versions), the reaction strength was revised to 3+ (Fig. 2).

DISCUSSION

In our study, the Eflexis analyzers could not determine the ABO group in 628 out of 19,466 cases. The Eflexis RR was 3.2% (including 1.1% for neonates up to 1 year of age), which was higher than that of 2.3% in 13,113 samples (including 1.2% for error messages) reported for the AutoVue system (Ortho-Clinical Diagnostics, Raritan, NJ, USA) [9]. However, in that study, samples from neonates up to 1 year of age were excluded, which may have affected the RR.

We did not use neonatal cards as these are more expensive than standard ABO RhD and reverse grouping cards. In total, 35.4% (222/628) of Eflexis rejected cases in our study were neonates. Anti-A and/or anti-B antibodies generally appear after the first 3–6 months of life, which explains why the AB group showed the lowest RR compared with the O group in neonates [15]. The RR depends on the proportion of neonates in each hospital and on whether a specific neonatal/cord blood ABO grouping card is available or whether only the standard ABO/Rh grouping card is used. Thus, laboratories should consider implementing a special interpretation process for neonates if only the standard ABO/Rh grouping cards are used. After further analysis, we discovered that 93.8% of the rejected reverse typing results in our laboratory were from neonates. Therefore, we applied a process to interpret ABO blood group results using only forward grouping and to disregard the reverse typing results for neonates of four months old or younger, as the reverse grouping results were not helpful for this group.

We found that the ABO group also affects the autoanalyzer RRs. Group A showed the highest RR, which can be explained by the highly frequent rejected results of reverse grouping because of the absence of expected reactions against B RBCs (Table 2). Likewise, weak plasma reactions on the AutoVue system causing ABO blood grouping discrepancies have been reported [9]. The number of A antigen sites ($0.81\text{--}1.17 \times 10^6$) per RBC is higher than that of B antigen sites ($0.61\text{--}0.83 \times 10^6$) [16]. This quantitative difference appears to be one of the reasons for the weaker reactions in reverse grouping using B RBCs than when using A RBCs.

Except for neonates, group O subjects showed the lowest RR (0.4%). Anti-A or anti-B in B or A subjects are mostly IgM in adults, but IgG anti-A and anti-B is far more common in O subjects [17]. We previously reported that ABO antibody titers were

higher in TUBE than in the gel card method using ID-Card (Dia-Med AG), except for the anti-A antibody in group O at the anti-globulin phase [18]. This result supports the concept that gel card reactions are weaker in IgM anti-A and anti-B of groups B and A than in IgG anti-A and anti-B of group O. However, the proportion of ABO groups varies among populations [19]. While group A is the most common group in most areas, most South American Indians have group O, and B is the most frequent group in Bengalese [19]. Therefore, the RRs of the Eflexis analyzer may vary among populations.

Hemagglutination reactions for reverse grouping by TUBE were significantly stronger than those for the autoanalyzers (Table 3); this phenomenon has also been reported for the AutoVue system [9].

“DP” was observed only in forward grouping for cases with weakened expression of A or B antigen on RBCs associated with ABO subgroups, neonates, or a mixture of A or B group RBCs that can be observed following ABO-mismatched transfusion or PBSCT, as reported in other studies [5, 9]. Therefore, when the Eflexis analyzer reports a “DP” flag, a weakened ABO group or mixture of different ABO groups should be considered. The “trace” and “?” flags are considered weak reactions in reverse grouping.

Depending on the co-inherited ABO allele, the *cis*-AB phenotype typically is A₂B₃, A₂B, or A₁B₃ [20]. Most *cis*-AB phenotypes have a weaker B phenotype than the normal B phenotype. Mistyping of *cis*-AB as normal AB was reported in 87.5% of cases when using the QWALYS-3 analyzer (Diagast, Loos, France) and 70.0% when using the Galileo NEO analyzer (Immucor Gamma, Norcross, GA, USA) [13]. The authors recommended the manual tile method as a simple supplemental test for the detection of the *cis*-AB phenotype, especially in countries with a relatively high prevalence of *cis*-AB phenotypes [13]. To prevent hemolytic transfusion reactions, O RBCs (or A RBCs when anti-A is not detectable in the plasma) and AB plasma or platelets are recommended for transfusion of patients with the *cis*-AB phenotype in Korea [21]. In our study, Aw or Bw other than *cis*-AB could be easily detected by the presence of “DP” flags and a weakened reaction in forward grouping. This can explain why there has been no problem with using the Eflexis prior to the current software upgrade as shown in case A in Fig. 2 in other countries. Since the software upgrade, the A₂B_w phenotype with 4+ strength in forward grouping using B RBCs has not yet been detected in our laboratory.

We also used arbitrary extra exclusion criteria that prevented automatic interfacing of the results. These criteria were designed to detect A₂B_w with weak reactions of anti-B antibodies. The au-

toanalyzers showed stronger reactions in forward grouping using anti-B anti-sera than TUBE for six cases of A₂B_w (Table 4). Our laboratory has a policy to repeat ABO blood grouping by TUBE following analysis on the autoanalyzer for all patients with an initial ABO blood grouping test. This policy and our extra exclusion criteria for forward grouping were helpful in detecting serologically suspicious ABO subgroups and preventing automatic interfacing of results without any flags from the autoanalyzers (case Nos. 8–13, Table 4). Although our laboratory policy for new patients is effective in identifying serologically suspicious ABO subgroups, additional criteria for forward grouping may help detect suspicious ABO subgroups that had previously been missed and grouped as a normal ABO group. Reviewing cases using the extra criteria regarding the hemagglutination reaction strength $\leq 3+$ for forward grouping, at least when using anti-B anti-sera, appears to be helpful in detecting ABO subgroups in regions where the *cis*-AB group occurs frequently [5]. However, the extra criterion of $\leq 1+$ hemagglutination reaction strength for reverse grouping did not help detect any extra plasma reactions of anti-B antibodies in four A₂B_w cases. Therefore, the extra criteria for reverse grouping other than those for forward grouping using Eflexis may be unproductive in increasing the probability of detecting A₂B_w, especially as it increases not only the laboratory workload but also the risk of clerical errors.

In summary, failed Eflexis analyzer ABO grouping was mainly due to various flags and discrepant ABO blood grouping with weak reactions in reverse grouping. This caused higher RRs in neonates and group A because our laboratory does not use neonatal cards; group O showed the lowest RR. As the *cis*-AB phenotype was overlooked by Eflexis, ABO grouping using TUBE in new patients and/or the implementation of extra rejection criteria in cases with $\leq 3+$ reaction grades for forward grouping were helpful in detecting the *cis*-AB phenotype. Following changes to the Eflexis software for scoring and the interpretation of weak hemagglutination reactions, there was a significant improvement in the detection of serologically suspicious ABO subgroups. The results of this study will help other laboratories to efficiently use Eflexis, especially in *cis*-AB phenotype-prevalent regions.

ACKNOWLEDGMENTS

We thank David Roxby for reviewing this manuscript.

AUTHOR CONTRIBUTIONS

Lim YA designed and supervised the study. Lim YA and Park SJ

wrote the manuscript. Lim YA, Park SJ, and Cho HS contributed to the data analysis and interpretation. All authors read and approved the final manuscript.

CONFLICTS OF INTEREST

No potential conflicts of interest relevant to this article were reported.

RESEARCH FUNDING

This study was supported by a Unionlab Research Fund in 2020.

ORCID

Young Ae Lim <https://orcid.org/0000-0002-8472-289X>
Seo-Jin Park <https://orcid.org/0000-0003-2768-1080>
Hyun Soo Cho <https://orcid.org/0000-0003-1086-2539>

REFERENCES

1. Elliott C, Carrascosa T, Souchet JL, Smith D, Malaxetxebarria I, You J, et al. Multicentre evaluation of Erytra Eflexis®, a benchtop fully automated analyser with a compact design for routine use in blood transfusion laboratory. *Transfus Med* 2019;29:401-7.
2. Roback JD, Denomme GA, Billingsley K, Bensing K, Parsons JC, McDonough WC. Performance and reliability of a benchtop automated instrument for transfusion testing: a comparative multicenter clinical study in the US population. *Transfusion* 2019;59:3511-8.
3. Chang C, Brown M, Davies L, Pointon L, Brown R, Barker D. Evaluation of Erytra® fully automated analyser for routine use in transfusion laboratory. *Transfus Med* 2014;24:33-8.
4. Taylor J, Hyare J, Stelfox P, Williams M, Lees R, Maley M. Multi-centre evaluation of pre-transfusional routine tests using 8-column format gel cards (DG Gel®). *Transfus Med* 2011;21:90-8.
5. Shin JW, Shin WY, Lee DL. Comparison of ABO blood group typing between automated blood bank analyzer IH-500 and manual method. *Korean J Blood Transfus* 2017;28:126-33.
6. Shin SY, Kwon KC, Koo SH, Park JW, Ko CS, Song JH, et al. Evaluation of two automated instruments for pre-transfusion testing: AutoVue Innova and Techno TwinStation. *Korean J Lab Med* 2008;28:214-20.
7. Cheng YW and Wilkinson JM. An experience of the introduction of a blood bank automation system (Ortho AutoVue Innova) in a regional acute hospital. *Transfus Apher Sci* 2015;53:58-63.
8. Morelatti F, Revelli N, Maffei LM, Poretti M, Santoro C, Parravicini A, et al. Evaluation of a new automated instrument for pretransfusion testing. *Transfusion* 1998;38:959-65.
9. Jo SY, Lee JM, Kim HL, Sin KH, Lee HJ, Chang CL, et al. Comparative analysis of clinical samples showing weak serum reaction on AutoVue system causing ABO blood typing discrepancies. *Ann Lab Med* 2017; 37:117-23.
10. Yamaguchi H, Okubo Y, Hazama F. An A2B3 phenotype blood showing atypical mode of inheritance. *Proc Jpn Acad* 1965;41:316-20.
11. Cho D, Kim SH, Jeon MJ, Choi KL, Kee SJ, Shin MG, et al. The serological and genetic basis of the cis-AB blood group in Korea. *Vox Sang* 2004; 87:41-3.
12. Han KS and Kang SH. Molecular genetic analysis of cis-AB blood group in Koreans. *Korean J Blood Transfus* 1997;8:85-91.
13. Chun S, Ryu MR, Cha SY, Seo JY, Cho D. ABO mistyping of cis-AB blood group by the automated microplate technique. *Transfus Med Hemother* 2018;45:5-10.
14. Fung MK, Eder AF, et al. (eds.) *Technical Manual*. 19th ed. Bethesda: American Association of Blood Banks, 2017: Methods 1-9.
15. Abo H and Le P. 1PK, GLOB, I and FORS blood group systems. In: Klein HG and Anstee DJ, eds. *Mollison's Blood Transfusion in Clinical Medicine*. 12th ed. Oxford: Wiley Blackwell, 2014:130.
16. Economidou J, Hughes-Jones NC, Gardner B. Quantitative measurements concerning A and B antigen sites. *Vox Sang* 1967;12:321-8.
17. Rawson AJ and Abelson NM. Studies of blood groups antibodies. III. Observations on the physicochemical properties of isohemagglutinins and isohemolysins. *J Immunol* 1960;85:636-9.
18. Kang SJ, Lim YA, Baik SY. Comparison of ABO antibody titers on the basis of the antibody detection method used. *Ann Lab Med* 2014;34: 300-6.
19. Abo H and Le P. 1PK, GLOB, I and FORS blood group systems. In: Klein HG and Anstee DJ, eds. *Mollison's Blood Transfusion in Clinical Medicine*. 12th ed. Oxford: Wiley Blackwell, 2014:119.
20. Chun S, Choi S, Yu H, Cho D. Cis AB, the blood group of many faces, is a conundrum to the novice eye. *Ann Lab Med* 2019;39:115-20.
21. Korea Centers for Disease Control & Prevention. *Korean National Transfusion Guideline*. 4th ed. Osong: Korea Centers for Disease Control & Prevention, 2016:15-6.